## Contact info

Aidian Oy aidian@aidian.eu / +358 10 309 3000 For orders: orders@aidian.eu

www.aidian.eu www.quikread.com www.tackleamr.com

QuikRead go CRP	Cat. no.
QuikRead go CRP 50 tests including capillaries and plungers	135171
QuikRead go CRP Control	153764
QuikRead go CRP Control High	153763

QuikRead go CRP+Hb	Cat. no.
QuikRead go CRP+Hb 50 tests including capillaries and plungers	140068
QuikRead go CRP Control	153764
QuikRead go CRP Control High	153763
QuikRead go Hb Control	141154

QuikRead go Strep A	Cat. no.
QuikRead go Strep A 50 tests including controls	135883

QuikRead go iFOBT	Cat. no.
QuikRead go iFOBT 50 tests	151051
QuikRead iFOB Sampling set 50 devices	151052
QuikRead go iFOBT Control Quantitative	152390
QuikRead FOB Positive Control	06027

QuikRead go Instrument	Cat. no.
QuikRead go Instrument (FI, SE, NO, DK, CZ, SK, PL, HU)	135867
QuikRead go Instrument (other languages)	133893


QuikRead go wrCRP	Cat. no.
QuikRead go wrCRP, 50 tests including capillaries and plungers	146521
QuikRead go CRP Control Low	153765
QuikRead go CRP Control	153764
QuikRead go CRP Control High	153763

QuikRead go wrCRP+Hb	Cat. no.
QuikRead go wrCRP+Hb, 50 tests including capillaries and plungers	146255
QuikRead go CRP Control Low	153765
QuikRead go CRP Control	153764
QuikRead go CRP Control High	153763
QuikRead go Hb 10 μl Control	153656

QuikRead go easy CRP	Cat. no.
QuikRead go easy CRP, 50 tests with 10 µl Sample Collector	153287
QuikRead go CRP Control Low	153765
QuikRead go CRP Control	153764
QuikRead go CRP Control High	153763
QuikRead go Sample Collector 10µl	153211

Information about QuikRead go accessories: www.aidian.eu

# It's so easy to use!







Result


- 1. Crocker J et al. Implementation of Point-of-Care Testing in an Ambulatory Practice of an Academic Medical Center. Am J Clin Pathol 2014; 142:640-646.
- 3. Minnaard M et al. Analytical performance. agreement and user-friendliness of five creactive protein point-of-care tests. Scand J Clin Lab Invest 2013; 73:627-634.
- 4. Hernández-Bou S et al. Evaluation of the bedside QuikRead go® CRP test in the management of febrile infants at the emergency department.



www.aidian.eu

# Results you can trust.

QuikRead go® for point of care





# Save time and resources with a flexible point of care system.





Available analytes

CRP / CRP+Hb / Strep A / FOB



Point of care testing (POCT)
can significantly improve
clinical operations with cost
reductions through improved
practice efficiency. The
economic benefits of POCT
may be realized in both feefor-service and global payment
environments."1

Get support for diagnosis and treatment decision with laboratory quality results just within minutes. Small and fully portable instrument enables testing in various healthcare settings.<sup>2-4</sup>



## **Immediate**

- Results in minutes
- Minimum hands-on time required
- Results you can trust during patient consultation
- Supports you on the spot in making the correct diagnosis and treatment decisions<sup>2-4</sup>



## Reliable

- Instrument read results ensure objective result interpretation
- Reliable results comparable to laboratory methods<sup>2-4</sup>
- Robust instrument
- Close to 50 000 QuikRead and QuikRead go instrument sold worldwide

Have peace of mind with the high performing, robust, and easy to use test system. Exceptional connectivity features and automatic testing procedure ensures you better POC management and traceability.<sup>2-4</sup>



# Carefree

- No need for regular maintenance
- Plug and play
- · Ready to use reagents
- Fully automatic testing procedure
- Interactive user interface available in 22 languages
- Battery use ensures full portability



## Traceable

- Results with all testing details stored automatically to the instrument
- Patient ID and User ID for improved traceability
- Connectable to most hospital (HIS) and laboratory information (LIS) systems
- POCT1-A2 bidirectional connectivity with advanced features
- LIS01-A2 unidirectional connectivity