

# CHELATING DETERGENT FOR THERMAL BEDPAN-WASHERS

## Anios LB 400

### CONCENTRATED DETERGENT RINSE FOR TREATMENT IN THERMAL WASH BASINS

- Detergent alkaline product
- High dispersing power of calcium salts
- Can be used in any quality of water (tap water or treated water)
- Chelating power
- Without EDTA nor NTA
- Phosphate free



### PRODUCT CHARACTERISTICS

- Clear to yellow liquid
- pH of pure product\* at +20°C: about 13.3
- pH diluted at 0.1%\*: about 11.1
- pH diluted at 0.3%\*: about 11.5
- Density (pure product) at 20°C: 1.13 g/ml
- Conductivity at 0.1% and 0.3%: 300 - 781  $\mu\text{S}/\text{cm}$
- Viscosity (pure product) : 3.84 mPa·s

\*pH value measurements under operating conditions always by means of a pH meter with temperature compensation

### PACKAGING AND ORDERING INFORMATION

PRODUCT	PACKAGING UNIT	ORDER CODE
Anios LB 400	2 x 5 L canister	XXXXXX



**LABORATOIRES ANIOS**  
1, Rue de l'Espoir  
59260 Lezennes  
France

Distributed by:  
**ECOLAB DEUTSCHLAND GmbH**  
Ecolab-Allee 1  
40789 Monheim am Rhein  
Germany

### APPLICATION AREA

Strong sequestering and dispersing product for use in hard to very hard water. Prevents limescale deposits in the device and on the equipment. Compatible with basins, glass, stainless steel and plastic urine bottles.

### INSTRUCTIONS FOR USE

Concentrated solution for professional use. Dilution according to the bedpan washer program and water hardness (1 to 3 ml/litre).

### INGREDIENTS

Sequestering agents in an alkaline medium.

### STORAGE


Dangerous. Respect the precautions for use (drawn up according to the European rules in force regarding the classification and labelling of chemical products). Storage: between +5°C and +25°C.

### CERTIFICATES

On demand we are pleased to provide you with certificates.

### MEDICAL DEVICE USE SAFELY

Please read label and product information.



FOR PROFESSIONAL USE ONLY.



www.ecolab.com